



WWW.FYICHINA.COM.CN

FANYUAN INSTRUMENT (HF) CO., LTD.
DOWELL SCIENCE AND TECHNOLOGY (HK) CO., LTD.
ADD: Tianyuan Road 9#, High-tech Zone, Hefei City, China Postal Code: 230088
TEL: +86-551-65232508
FAX: +86-551-65232507
Email: sales@hkdownell.com
Website: http://www.fyichina.com.cn

GD350 GARMENT SAMPLE DYEING MACHINE





FANYUAN INSTRUMENT (HF) CO., LTD.
DOWELL SCIENCE AND TECHNOLOGY (HK) CO., LTD.
ADD: Tianyuan Road 9#, High-tech Zone, Hefei City, China Postal Code: 230088
TEL: +86-551-65232508
FAX: +86-551-65232507
Email: sales@hkdownell.com
Website: http://www.fyichina.com.cn

Model No.

GD350

Application:

GD350 garment sample dyeing is suitable for the dyeing of all kinds of garment, such as linen, cotton, rayon, etc.

Features:

- Structure of mirror stainless steel;
- Energy- saving type but multiple function
- Machine with lower liquor ratio
- Easy to operate and computer control
- reversible electric motor and reduction system
- Independent and high precision heating system
- Adjustable rotating-speed drum with low noise;

Key Specification:

Liquor ratio	1: 10—1: 20
Drum Size	∅ 350mm × 600mm
Water tank	Approx. 650 × 400 × 550mm
Rotating speed	0-40rpm, reversible rotation
Power of motor	250W
Temperature sensor	PT-100 thermotic couple
Working temperature	RT~99℃
Temperature accuracy	±1℃
Power supply	220V AC, 50Hz
Dimensions	1150×470×820mm
Net weight	80kg



FANYUAN INSTRUMENT (HF) CO., LTD.
DOWELL SCIENCE AND TECHNOLOGY (HK) CO., LTD.
ADD: Tianyuan Road 9#, High-tech Zone, Hefei City, China Postal Code: 230088
TEL: +86-551-65232508
FAX: +86-551-65232507
Email: sales@hkdownell.com
Website: http://www.fyichina.com.cn

Standard configuration:

No.	Item	Quantity
1	Main machine	1 set

Optional Accessories:

No.	Item